



MORPHEUS ND

The New Dimension of Anaesthesia



we are glad to introduce the latest innovation in Health Care Technology as a result of 44 years of hi-specialized R&D.

MORPHEUS ND, the smaller and lighter Anaesthesia workstation existing in the world with a main 15" touch screen display including all the advanced anaesthesia parameters, suitable for adults, children and neonatal patients.

The unit incorporates the HYBRID exclusive system with a silent low pressure turbine for the continuous working of the lung ventilator even in absence of compressed air and oxygen. In the event of turbine failure the unit can work with compressed gas as an additional feeding source.

A full programme of safety devices is also present as well as the USB output. Among them the exclusive PROTOLOCK patented system against the accidental exchange of O₂/N₂O gases during anaesthesia.

Morpheus ND can be mounted on a pendant unit in operating room or on a simple trolley.

The unit is 100% made in Italy, UE.